



time to relax



SKYLINER 25Y



< timetorelax

True relaxation begins not with a button – but with trust.

For over two decades, CASADA has been synonymous with refined wellness technology – offering more than function: a promise of comfort, craftsmanship and ease.

Blending innovation with timeless design, our products are made for those who seek lasting wellbeing, aesthetic harmony and everyday indulgence.

With a presence across the globe, CASADA remains devoted to one guiding principle: excellence in every detail.

Congratulations

With this choice, you've shown how much you value your wellbeing and the art of relaxation. To enjoy your chair for years to come, we kindly ask you to read the safety instructions carefully and follow them.

May your new favourite place bring you countless moments of calm, comfort, and quiet joy.

Table of Contents

Safety Instructions	04
Contraindications.....	05
Warnings on Use	06
Features Functions	07
Set-up and Assembly.....	08
Cleaning Maintenance Storage	14
Structure of the remote control.....	16
Operate massage chair	18
Automatic programs.....	21
Remote Control.....	25
Adjusting the Chair	36
Settings	38
Other functions	41
Memory.....	46
Troubleshooting	47
Specifications.....	47
EU–Declaration of Conformity	48
Warranty Conditions.....	48

Safety Instructions



Please read these instructions carefully before using the massage chair. Doing so will help ensure safe operation and extend the life of your device. Keep this manual in a safe place for future reference

- This massage chair complies with applicable technical standards and current safety regulations.
- The unit includes heated areas. Users with sensitive skin or reduced heat sensitivity should exercise caution when using the product. (ICE 60335-2-32)
- The massage chair is maintenance-free. Repairs must only be carried out by authorised service personnel.
- Improper use and unauthorised repairs are prohibited for safety reasons and will void the warranty.
- Take care not to insert fingers or objects between the massage rollers.
- Never touch the power plug with wet hands.
- Avoid exposing the device to water, direct sunlight or extreme temperatures.
- Do not use the massage chair in damp or humid environments (e.g. bathrooms).
- Ensure that ventilation openings are not blocked and that adequate airflow is maintained.
- Disconnect the massage chair from the mains power supply when not in use for an extended period.
- Do not use damaged plugs, cables or sockets.
- Damaged power cords or plugs must be replaced by the manufacturer, an authorised service centre or qualified personnel.
- In case of malfunction, disconnect the device from the power supply immediately.
- Using the device for purposes other than those intended may result in damage and void any liability.
- For your safety and comfort, do not exceed a daily usage time of 30 minutes.
- Avoid using sharp or pointed objects on or near the device.
- Use only suitable sockets and ensure the plug is fully inserted.
- Always disconnect the device from the power supply before cleaning or carrying out any maintenance.
- Children under the age of 14 may only clean or operate the device under adult supervision.
- In the event of a power failure, switch the unit to "OFF" and unplug it.
- Do not operate the massage chair if the cable or plug is wet or damaged.
- Children under 14 and individuals with limited physical, sensory or mental capabilities must only use the device under supervision.
- Do not use the massage chair immediately after eating or while under the influence of alcohol.
- Use the device only as described in this manual.
- Only use original accessories provided by the manufacturer.
- Do not stand, sit or jump on the footrest.
- Do not sit on the armrests or backrest.
- Keep the massage chair clean and ensure that no foreign objects fall into any gaps or openings.
- Never pull directly on the power cord.
- When adjusting the footrest or backrest, make sure there are no objects or people underneath – especially children.
- This product is intended for private indoor use only.
- If the power cord is damaged, it must be replaced with an original replacement cord provided by the manufacturer or an authorised service provider.

Contraindications and Important Usage Information

If you are unsure whether the massage chair is suitable for you, please consult your doctor before using the device.

Do not use the chair if you are pregnant or have any of the following conditions in the massage area:

- recent injuries
- thrombotic disorders
- acute inflammation or swelling
- existing cancer diagnoses
- If you are undergoing medical treatment or have any health concerns, prior consultation with your doctor is strongly recommended.
-

Contraindications

Usage Restrictions

Individuals with physical, sensory, or mental impairments, as well as those lacking experience or sufficient understanding, should only use the massage chair under the supervision and guidance of a responsible person.

Medical advice must be obtained before use in the following cases:

- individuals who are bedridden
- those undergoing ongoing medical treatment
- known back conditions or spinal disorders
- use of a pacemaker

Usage Instructions

- Stop using the chair immediately if you experience any pain.
- If an injury occurs during use, seek medical attention without delay.
- Do not massage reddened, swollen, or inflamed areas of skin.
- This massage chair is intended solely for private wellness use. It is not a medical device and does not replace professional treatment.
- The recommended maximum daily usage time is 30 minutes. Do not massage the same area for more than 15 minutes to avoid overstimulation or discomfort.
- This product is not suitable for children under the age of 14. Children may only use the chair under supervision and must not treat it as a toy.

Contraindications for the Use of BRAINTRONICS®

The use of the braintronics® programme is not recommended in the following cases:

- Psychotic disorders (e.g. schizophrenia, bipolar disorder, endogenous depression)
- Epilepsy or other seizure-related conditions
- Heart conditions or recent heart attack/stroke
- Disorders of the central nervous system
- Thrombosis
- Attention deficit disorders (inattentive type)
- Personality disorders
- Intellectual disability
- Addiction-related conditions (e.g. alcohol, drug, or medication dependency)
- Regular use of psychotropic medication
- Pregnancy

Important Notice:


For ethical and legal reasons, braintronics® must not be used by children or adolescents without the explicit consent of a legal guardian.

Individuals who may experience inner or religious conflict because of meditative states are advised to refrain from using the programme.

Warnings on Use


NOTE

For your own safety and to avoid damaging the massage chair, do not exert pressure on the footrest when extended.





NOTE

Before each massage, make sure that the massage rollers are not in the seat area.



Features | Functions

1. Pillow
2. 3D digital speakers
3. Shoulder airbag
4. Backrest cushion
5. Lumbar airbag
6. Ion air blower
7. Quick launch buttons
8. USB charging connector
9. Arm airbag
10. Storage space
11. Seat cushion
12. Calf massage
13. Armrest
14. LED lamps
15. Foot massage



Set-up and Assembly

1. Assembling the armrest

- From the arm rest, first remove the connectors and the air hose which match those on the frame of the massage chair (Fig. 1).
- Now link each connection with its counterpart (colour-coded). (Fig. 2).




Fig. 1




Fig. 2

- Insert the connecting parts on the underside of the frame firmly into the armrest connection. Correctly position the compressed air line. (Fig. 3).
- After connecting with the appropriate connection lines, carefully insert the middle of the steel pipe behind the steel frame manually into the lower part of the frame, making sure that the compressed air line is not blocked or damaged. (Fig. 4).




Fig. 3




Fig. 4

- Now hang the bolt on the rear side panel (inner side) into the guide on the massage chair. (Fig. 5)
- Lock the front attachment of the seat in position and use slight pressure to hook this into place. (Fig. 6)




Fig. 5




Fig. 6

- Use 2x (M5X10) screws to secure the front attachment. (Fig. 7).
- Use 5x (M4x10) screws to secure the rear side of the side panel. (Fig. 8).




Fig. 7




Fig. 8

- Narrow side first, hook the covers into the side panel and then clip the mask onto the side panel. (Fig. 9).
- From the inner side of the shoulder part, use one ST4,2x16 screw to secure the cover. When assembling, the assembly instructions for the two side panels are identical. (Fig. 10)




Fig. 9




Fig. 10

Set-up and Assembly

2. Assembling the footrest

- Position the foot rest in front of the seat and connect the air hose and plug to the counterpart on the seat (Figs. 1 + 2).
- Then push any protruding air hose and cable back into the seat.




Fig. 1




Fig. 2

- Slide the inner and outer sleeves onto the mounting bracket on the footrest. Check that the head end is facing the footrest (Fig. 3).
- On the inner side of the seat, open the hinges on the right and left (Fig. 4).




Fig. 3




Fig. 4

- Hang the footrest with the sleeves into the hinges (Fig. 5)
- Close the hinges and secure them using an M5 screw. (Fig. 6)




Fig. 5




Fig. 6

- Then place the plastic covers over the hinges. (Fig. 7)




Fig. 7

3. Assembling the touchscreen bracket

- Insert the holder into the side panel and secure it using 3x M4 screws. (Figs.1+2).




Fig. 1




Fig. 2

Set-up and Assembly

- Click the plastic foot of the bracket onto the installation surface of the massage chair. (Figs. 3+4)




Fig. 3





Fig. 4

NOTE: The connection for the remote control is located on the inside of the holding under the seat cushion. Insert the cable and secure it.

4. Using the pillow and backrest cushion


- The kneading intensity of the massage on the neck and shoulder can be alleviated by using the pillow. We recommend using the cushion.
- The backrest cushion is connected to the backrest by a zipper (1), and the pillow is connected to the backrest cushion by a Velcro strip (2).



Set-up and Assembly

5. Floor protection




- To avoid damaging the floor, place the massage chair on a carpet or any other stable surface.



CAUTION: Make sure that you do not move the chair across wooden floors, uneven floors or in confined spaces.

6. Switching off, finishing the massage


- By pressing the power switch (1) during a massage, all massage functions are switched off immediately. The backrest and calf unit are reset, the massage time expires and all massage functions are switched off immediately.
- The power supply throughout the entire massage chair is switched off. (2) (3)



Set-up and Assembly

7. Repositioning the massage chair

- Make sure no cables are lying across the floor. Tilt the backrest somewhat backwards (with the centre of gravity on the wheels). Use your hands to push the chair forwards or backwards, then return it to its original position slowly and gently.



NOTE: Before moving the chair, you must first switch it off and then disconnect the power cord and the connecting cable.

NOTE: Please do not move the chair alone. The second person can lift the chair at the front to position it onto the rollers.

Cleaning | Maintenance | Storage

A | CLEANING THE SYNTHETIC LEATHER

Wipe the synthetic leather with a soft, dry cloth. Please do not clean with products which contain chemicals, such as thinners, benzine and alcohol. Make sure that the cleaning agents do not damage the synthetic leather. Please read the instructions on the cleaning agent very carefully. If the synthetic leather is very dirty, clean it as follows:

1. Moisten a soft cloth with hot water or a commercially-available neutral detergent and wring it out.
 2. Check that the cloth has been well wrung out.
 3. Clean the surface of the synthetic leather with the cloth.
 4. While cleaning, make sure that you rinse the cloth several times.
 5. Leave the synthetic leather to dry naturally. Do not use a hairdryer.
- The surface of the synthetic leather may fade or discolour if an unsuitable cleaning agent is used. Before using any cleaning solution on the entire surface, test it first on a small section.
 - Cleaning the synthetic leather with alcohol, thinner or petrol-containing agents may lead to discoloration or even remove the paint.


B | CLEANING THE PLASTIC PARTS

1. Put a soft cloth in warm water with a neutral detergent, wring out well and wipe off the plastic parts.
2. Do not use chemicals, such as thinners, benzine and alcohol to clean the plastic parts.
3. Wipe with a microfibre or soft cloth wrung out in clean water.
4. Allow the plastic parts to air dry.

C | STORAGE

- If the device is not used for a longer period of time, the power cable should be wound up and the device should be stored in a dry environment.
- Never store the device at high temperatures or in proximity to open fire. Avoid prolonged exposure to sunlight.
- When cleaning, only use commercially-available products for synthetic leather. Never clean it with thinner, benzene, or alcohol.
- If the device is on an uneven surface, do not roll or push it; lift it up in order to move it.
- Please take breaks after each use and do not use the device a number of times without pauses.

Structure of the remote control



1. On-/ Off-Button
2. Pause or continue the massage
3. Activates or deactivates the belt heat
4. Activates or deactivates the O2 ion function
5. Activates or deactivates the foot rollers
6. Activates or deactivates calf knead function
7. Adjusting the massage intensity
8. Turns the massage chair on or off
9. Adjusting the air pressure for the lower body
(Only possible in manual mode)
11. Adjusting the shoulder position


Structure of the remote control




12. Display of the current status of the massage chair
13. Selection of automatic massage programs
14. Selection of manually adjustable functions
15. Adjusting the position of the seat and the angle of the sitting posture
16. Settings
17. Selection key for the massage techniques
18. Increase or decrease the air pressure intensity on the upper body
19. Status display
20. Selection of the compression massage area
21. Selection of the length of the massage program
22. Adjusting the massage speed
23. Adjusting the 4D intensity
24. Custom memory

Operate massage chair

1. Press the On/Off button on the side of the remote control  or long press the On/Off button on the armrest for the system to start.




2. After starting the system, select the massage intensity by pressing the "Gentle"  "Soothing"  "Strong"  buttons on the remote control. This switches on the programme, and the massage chair changes to the lying position. If no selection is made within five minutes, the chair will automatically switch off.




3. For the massage programmes, select the automatic programme (see 7.3) and the massage chair will be moved to the zero gravity position. Wait until the setting process is finished. Press the Skip button , to skip a step and go directly to the next level.

Operate massage chair


4. The chair now scans your body shape. Wait until your body shape check has been fully scanned. You can also press the Skip button  or skip the test and start the massage directly.



After scanning your body shape, the remote control will give an acoustic signal. The signal alerts you that you can now fine-tune the shoulder position. When the massage roller is in the right position for you, press the OK button , to start the massage. If the position is not correct, use the Up button  or the Down button , to manually correct the position of the massage and press the OK button  to confirm your selection and start the massage. If no buttons are pushed within 20 seconds, the rollers will move to the starting position and the massage will start.



Operate massage chair



Pause function

(Please use this function only after switching on the chair [see page 77].)


After switching on the massage chair, press the Pause button  to pause the massage process. Press the Pause button  again to resume the programme. Please note: In pause mode, all functions except shutdown are disabled.

Automatic programs

(Please use this function only after switching on the chair [see page 77].)

Selecting the Automatic programs

1. Select the Auto button  in the left menu bar.
2. Select the desired automatic massage programme and click to start the massage.



You can choose from 20 automatic programmes.

Automatic programs

AI Intelligent - Programs	
Function	Description
Total Relax	The upper mechanical core is used to massage the neck, shoulders and Fengchi acupoint by kneading. This is followed by tapping and shiatsu to massage the whole body. The lower mechanical core is used to massage the waist and hips by kneading and tapping. The two mechanical cores coordinate to help relax the muscles of the entire body and relieve fatigue.
Stress Release	Massages the neck, shoulders and Dazhui acupuncture point mainly by kneading to relieve systemic pressure.
Full Airbag	Applies compression massage to the whole body, together with the sole roller and calf kneading massage.
Chinese Massage	Using the massage technique of Traditional Chinese Medicine, 4D massage is applied to the back and waist. Thus, the blood circulation of the back and waist is stimulated, the meridians are soothed and blockages in the back and waist are released by the massage.
Thai Massage	This massage dominates with the help of compression massage on legs, hands, shoulders and hips. During the movement, the 3D massage hand holds with pressure against the waist to cause the Thai stretching. Note: The stretching is divided into three levels. For a gentle stretch, select the lowest level. For a strong stretch, select the strongest level.

Automatic programs

Relaxing and activating - Programs	
Function	Description
Meridian Treatment	Massage of the shoulders, neck, back and waist using kneading, tapping and shiatsu techniques to treat the meridians.
Neck & Shoulder	Massages the neck and shoulders mainly with kneading and stretching techniques and the upper back with a medical massage. Jianjing acupuncture points are pressed to relieve cervical spondylosis and relax the muscles of the neck and shoulders.
Waist Care	Massages the waist by kneading, kneading&tapping and tapping. This is followed by a 4D massage at the waist.
Waist & buttock Shape	Massages the buttocks through kneading and a medical massage. This is followed by a deep massage for the buttocks and perineum through tapping and kneading. To relax the hip muscles, raise your hip.
Leg & foot Relieve	Massages the shoulders and waist mainly by kneading, tapping and kneading&tapping. Applies a soothing compression massage to the calves and feet.
Health Recharge	Deep massage of the whole body using kneading, tapping and shiatsu techniques to promote blood circulation throughout the body and prevent blood clots. It is especially for people who do not get much exercise. The massage intensity is moderately gentle.

Automatic programs

Exclusiv for Family - Programs	
Function	Description
Power	Massage the shoulders, neck, back and waist mainly by kneading, tapping and kneading&tapping. This relieves fatigue, promotes blood circulation and reduces stress through deep massage. The massage intensity is moderate, making it suitable for men.
Keep fit	Massages the shoulders, neck, waist and hips through kneading and shiatsu techniques. By massaging the whole body, it causes the body and mind to relax. Enjoy treatment worthy of a queen. This is particularly suitable for women, as the massage intensity is very gentle.
Senior	Applies a gentle massage to the whole body in addition to the sole roller and calf kneading massage.
Office Regimen	Massages the shoulders, neck and waist with the massage techniques of kneading, tapping and kneading&tapping. Deep massage is used to relieve shoulder and neck pain, relieve sore muscles in the lumbar region and maintain the lumbar spine. It is especially for office workers and for people who sit a lot.

Automatic programs

Overall Relaxation - Programs	
Function	Description
Energy Boost	Deep massage of the whole body using kneading, tapping and shiatsu techniques to awaken the vitality of the human body and stimulate energy.
Lunch Break	Massages the back and waist using kneading and tapping techniques while pressing the Shenyu acupoint. This massage aims to improve the quality of your lunch break so that you can start the afternoon with a new burst of energy.
Sweet Dream	Massages the shoulders, neck, back and waist through kneading and a medical massage technique. The Fengchi and Shenyu acupuncture points are pressed to improve sleep quality and promote sleep. It is especially for people who suffer from sleep disorders.
Brain Refresh	Massages the shoulders, neck, back and waist through kneading, 4D kneading and Shiatsu techniques. A high degree of physical and mental harmony is achieved through deep massage. It is especially for people who are mentally exhausted and overworked. The intensity of the massage is gentle.
Rejuvenate	Massages the back and waist, mainly with massage techniques such as tapping and kneading&knocking. This promotes blood circulation, eliminates local muscle pain, accelerates the elimination of lactic acid and regenerates human vitality.

Remote Control

Program selection

Select the  button on the left menu bar to enter the selection menu. You can choose from the following settings: Massagetechnique, Massage Area, 4D intensity, Position, Massage Width und Massage Speed.

- There are a total of 11 options in the massage technique menu: Stopp, Knead, Tap I, Tap II, Knead & Tap, Shiatsu I, Shiatsu II, Grap, 4D I, 4D II, 4D III, 4D IV. (Please note: To display more massage techniques, press the arrow buttons  or  on the remote control)




Function	Setting the massage rollers	
4D I 	Speed and strength	- Speed can be set at five levels. - Width can be set at five levels. - Intensity can be set at five levels.
4D II 	Speed and strength	
4D III 	Speed, width and strength	
4D IV 	Speed, width and strength	

Remote control

Function	Setting the Massage Rollers	
Stop		Press this button for massaging without a technique.
Knead		Speed and strength can be set.
Tap I		Speed, width and strength can be set.
Tap II		Speed, width and strength can be set.
Knead & Tap		Speed, width and strength can be set.
Shiatsu I		Speed, width and strength can be set.
Shiatsu II		Speed, width and strength can be set.
Grasp		Speed and strength can be set.

- Speed can be set at five levels.
- Width can be set at five levels.
- Intensity can be set at five levels.

Please note: You can also press the Massage Technique button  on the “Home” screen  where you can select the desired massage technique from the pop-up menu. (This function is a shortcut for changing massage techniques on the current “Home” screen of the remote control programmes.)



Remote control

2. The massage area menu offers you five areas as options: Point-by-Point, Partial, Full, Upper (back) and Lower (back). Select the desired area for the massage.



Button	Description
	Fixed point massage
	Local forwards and backwards massage in a partial area.
	Full forwards and backwards massage.
	Forwards and backwards massage of the upper back.
	Forwards and backwards massage of the lower back.

Please note: You can also click on the massage area button  on the “Home” screen  to select the corresponding massage area in the pop-up menu. (This function is a shortcut for switching the massage areas on the current “Home” screen of the remote control programme.)

Remote control



3. You can set the 4D intensity in the menu. Setting levels 1-5 are available, with 1 being the lowest and 5 the strongest. The higher the 4D intensity is set, the further the massage roller protrudes and the stronger the massage. Select the Plus **+** button or Minus button **-**, to adjust the 4D intensity.



Remote control

Level 0	Level 1	Level 2	Level 3	Level 4	Level 5

Button	Description
	Press this button to increase the 4D intensity in levels 1 to 5.
	Press this key to decrease the 4D intensity in levels 1 to 5.

Please note: You can also set the 4D intensity with the Plus button **4D+** or the Minus button **4D-** in the "Home" screen. (This function is a shortcut for setting the 4D intensity on the current "Home" screen of the remote control programme.)



Level 0	Level 1	Level 2	Level 3	Level 4	Level 5

Button	Description
	Press this button to increase the 4D intensity in levels 1 to 5.
	Press this button to decrease the 4D intensity in levels 1 to 5.

Remote control

4. To set the massage robot, select the arrow pointing upwards ▲ or the downward pointing arrow ▼, to fine-tune the massage hand up or down.



5. Select the message widths as shown in the illustration below.



- (1) The width can only be set when the user activates the remote control mode when selecting the massage technique or massage area.
- (2) The massage width can be set from level 1 to 5, where 1 is the narrowest and 5 is the widest setting. Please select one of the following massage techniques: Tap I, Tap II, Shiatsu I, Shiatsu II, 4D III und 4D IV.

Remote control

Level 0	Level 1	Level 2	Level 3	Level 4	Level 5
---------	---------	---------	---------	---------	---------

Button	Description
+	Press this button to increase the massage width in steps 1 to 5.
-	Press this button to decrease the massage width in steps 1 to 5.

6. Select the message speed as shown in the illustration below.



- (1) The speed can only be set if no automatic massage programme is running at the same time.
- (2) The massage speed can be set from level 1 to 5, where 1 is the slowest and 5 is the fastest setting. The speed cannot be adjusted if the massage process has been stopped.

Level 0	Level 1	Level 2	Level 3	Level 4	Level 5
---------	---------	---------	---------	---------	---------

Button	Description
+	Press this button to increase the massage speed in levels 1 to 5.
-	Press this key to decrease the massage speed in levels 1 to 5.

Please note: You can also use the massage speed + **SPEED+** or massage speed - **SPEED-** on the "Home" screen 🖱️ to set the massage speed.

Remote control




Level 0	Level 1	Level 2	Level 3	Level 4	Level 5
---------	---------	---------	---------	---------	---------

Button	Description
SPEED+	Press this button to increase the massage speed in steps 1 to 5.
SPEED-	Press this button to decrease the massage speed in steps 1 to 5.

Remote control

7. Go to the Manual menu to select the Air mode as shown in the following figure. The air mode includes full back, upper back or lower back. There is also a hand roller massage function on the inside of the arm that automatically turns on when the full body air massage is selected.



Button	Description
	Press this button to activate or deactivate the <i>complete back</i> air massage.
	Press this button to activate or deactivate the <i>upper back</i> air massage.
	Press this button to activate or deactivate the <i>lower back</i> air massage.

Air massage display		
Complete body air pressure	Upper body air pressure	Lower body air pressure

Remote control

Please note: You can also press the air massage button  on the "Home" screen  and select the desired air massage area from the pop-up menu. (This function is a quick key to the air massage on the current "Home" page of the remote control).



8. To set the air intensity, the air massage must be switched on. There are a total of 5 levels, with 1 being the lowest and 5 the highest level.



Remote control

					
Level 0	Level 1	Level 2	Level 3	Level 4	Level 5

Button	Description
	Press this key to increase the compression intensity between 1 to 5.
	Press this key to decrease the compression intensity between 1 to 5.

Please note: You can select the Upper Body Air Massage + button **Upper Air +** or Lower Body Air Massage - **Upper Air -** on the "Home" Screen  or set the upper body compression intensity; the settings can also be set in the same manner with the Lower Body **Lower Air +** **Lower Air -**. There are 5 levels of compression intensity.



					
Level 0	Level 1	Level 2	Level 3	Level 4	Level 5

Button	Description
Lower Air + Upper Air +	Press this button to increase the compression intensity between 1 and 5.
Lower Air - Upper Air -	Press this button to decrease the compression intensity between 1 and 5.

Adjusting the Chair

Click on the Angle button  in the menu bar on the left side of the screen to enter the settings menu where you can adjust the angle of the chair. Adjustments include extending the legrest, extending the footrest, raising the legrest, raising the joint, raising the backrest, zero gravity I and II and other setting functions.



1. Setting zero gravity

Button	Description
Zero Gravity I 	Zero Gravity I button: automatically setting the legrest and backrest position to Zero Gravity I position.
Zero Gravity II 	Zero Gravity II button: automatically setting the legrest and backrest position to Zero Gravity II position.

2. Lengthening or shortening the footrest

Button	Description
Shortening footrest 	Button for shortening the footrest: To manually set the footrest, press and hold this button to slowly move the footrest upwards and shorten the footrest. Release the button when you are finished with the manual setting process.
Lengthening footrest 	Button for lengthening the footrest: To manually set the footrest, press and hold this button to slowly move the footrest downwards and lengthen the footrest. Release the button when you are finished with the manual setting process.

Adjusting the Chair

3. Lengthening or shortening the legrest

Button	Description
Shortening legrest 	Button for shortening the legrest: Shorten the legrest using the remote control. Press and hold this button to slowly move the legrest upwards. Release the button to stop.
Lengthening legrest 	Button for lengthening the legrest: Lengthen the legrest using the remote control. Press and hold this button to slowly move the legrest downwards. Release the button to stop.

4. Raising and lowering the legrest

Button	Description
Raising legrest 	Raising the legrest: Press and hold this button to slowly raise the legrest. Release the button to stop the movement. When the legrest has been fully raised, the footrest is automatically adjusted by the chair.
Lowering legrest 	Lowering the legrest: Press and hold this button to slowly lower the legrest. Release the button to stop the movement. When the legrest has been fully lowered, the footrest is automatically adjusted by the chair.

5. Setting the angle

Button	Description
Raising chair 	Raising the chair: Press and hold this button to slowly raise the chair angle and release the button to stop the movement of the backrest upwards movement or the legrest downwards. When the desired angle has been reached, the footrest is automatically adjusted.
Lowering chair 	Lowering the chair: Press and hold this button to slowly lower the joint while the legrest is slowly raised. Release the button to stop the movement of the backrest downwards or the legrest upwards. When the desired angle has been reached, the footrest is automatically adjusted.

6. Setting the back angle

Button	Description
Raising backrest 	Button for raising the backrest: Press and hold this button to slowly raise the backrest. Releasing the button again stops the function.
Lowering backrest 	Button for lowering the backrest: Press and hold this button to slowly lower the backrest. Releasing the button again stops the function.

Settings


Tap the  button in the menu bar on the left side of the screen to access the settings. This is where the model, software version, serial number and firmware version are displayed. The interface also includes the language settings.

1. The Current Unit information shows the model, software version and the serial number of the product and firmware.
2. Bluetooth (playback via Bluetooth speakers)
Once the sound source (such as a mobile phone, MID tablet, etc.) has been paired with the Bluetooth function of the massage chair, the music will be transmitted wirelessly to the playback system of the massage chair.
 - a. Access the main menu of your smart device (e.g., mobile phone). Open Settings and select the Bluetooth option. Activate the Bluetooth and search for devices in your vicinity.
 - b. The smart device will search for the corresponding Bluetooth name - Skyliner II. Connect your device to the chair.
 - c. Now start the desired music on your device.

Settings

Reset


Select whether the chair should automatically return to its original position after the massage or not.



Settings

5. Language selection

The standard language is English. The following languages can be selected: English, German, Hungarian, Polish, French, Russian and Spanish.



Other functions

On the display under "Home" , select the "Belt heater" button , to switch the heating function of the wrap (graphene heating) on or off. After switching on the heating function, heat will be emitted for 3 minutes.



Other functions

On the display under “Home” , select the “O2 Ion” button , to switch the O2 ion heating function on or off.



On the “Home” screen , click on the “Roller” button , to switch the foot massage function on or off with three adjustable speeds.



Button			Description
			Press this button to activate or deactivate the sole roller function. Press the cyclically programmed button multiple times to switch between speed levels 1, 2, 3 and off.
Level 1	Level 2	Level 3	

Other functions

On the display under “Home” , select the “Calf” button , to switch the calf massage function on or off. There is a friction massage function for the inside of the calf and an air cushion massage function for the outside, which can relieve stiffness in the calf muscles or tension caused by standing for long periods of time.



Other function

On the display under “Home” , select the “Shoulder” button  and press the up and down arrow buttons in the pop-up menu to adjust the shoulder position to a total of 11 levels.



Other functions

On the “Home” screen display , select the “Massage Time Setting” button . Press the plus and minus arrow button in the pop-up menu to lengthen or shorten the massage time. The time is always shortened or lengthened in 5 minute steps. The longest running time is 40 minutes.



On the “Home” screen display , select the “Soft” , „Soothing”  or „Strong”  button to obtain the desired massage intensity.



Memory

Click on the “Memory” button  in the menu bar on the left side of the screen to access the massage function memories. There are two categories: saving and recall of the massage function. The massage can be saved under “Save I”, “Save II” or “Save III”. You can call up the saved massage under “Use I”, “Use II” and “Use III”.

Button	Description
	Saving: Save the currently running massage programme, including automatic programme, massage technique, air mode and intensity, position and angle of the leg and backrest unit and other massage functions.
	Recall: Recall the previously saved massage program, including all settings.



Troubleshooting

FAULTS	EXPLANATION
• The noise of the engine is audible during use. • The remote control does not work properly. • The device stops operating. • The braintronics® programmes are not starting.	• This is a normal operating noise. • Ensure that the mains plug is properly inserted into the socket and that the power switch is switched on. Please check the cable of the remote control. • The chair switches itself off automatically once the pre-set runtime has been reached or after excessive usage. This ensures the chair has a long operational life. • Make sure that you have installed the braintronics app and activated your account. Also check there is a Bluetooth® connection between the smartphone and the massage chair.

Specifications

Dimensions:	upright: 80x148x120 cm reclined: 80x180x105 cm
Weight:	141 kg
Voltage:	220-240V ~ 50/60Hz
Nominal power:	220 W
Runtime:	20 minutes
braintronics:	21 minutes
Certificates:	CE X EAC

EU – Declaration of conformity

Complies with European Regulatory Standards EMC 2014/30/EU and LVD 2014/35/EU.

Warranty conditions

The retailer provides the statutory warranty period for the product described on the reverse side. This warranty period begins on the date of purchase, which must be verified with a valid proof of purchase. Within Germany, the retailer will repair or replace products found to be defective during the warranty period at no cost to the customer, provided the fault is not due to improper use or external influence. Components subject to normal wear and tear, such as fabric covers, are excluded from the warranty.

To claim warranty service, the defective product must be returned along with proof of purchase before the end of the warranty period. Warranty claims are void if the fault was caused by external factors, improper handling, or by repairs or modifications not carried out by the manufacturer or an authorised service partner.

The warranty is limited to the repair or replacement of the product. Any further claims, in particular compensation for damages, are excluded within the limits of applicable law. If the manufacturer's inspection reveals that the product fault is not covered under the warranty or that the warranty period has expired, the customer shall bear the cost of inspection and any repair.

Legal Disclaimer: Use and Health-Related Information

All information provided in connection with this product and its applications does not constitute medical advice or therapeutic claims. The diagnosis and treatment of illnesses and other physical disorders require consultation with a licensed medical professional, therapist, or alternative health practitioner. The content is intended solely for informational purposes and must not be used as a substitute for medical advice or treatment. Each user is responsible for determining the suitability of the application in their individual case and is advised to consult a specialist if in doubt.

All use of the product is at the user's own risk. We expressly disclaim any healing claims or health-related promises.

Casada International GmbH
Obermeiers Feld 3
D – 33104 Paderborn Germany
info@casada.com
www.casada.com



Copyright Images and text are subject to the copyright of CASADA International GmbH and may not be used without express consent.

Copyright © 2021 CASADA International GmbH.
All rights reserved